



2441 Cades Way, Suite 200
Vista, CA 92081
Phone (888) 310-7374
Fax (866) 787-7597
www.uc-scanning.com

REMOTE LINK INSTALL INSTRUCTIONS:

1. Insert CD into the CD-ROM drive
2. Wait for the auto-start setup window to pop up in the Windows Operating System. If the auto-start window does not appear, follow instructions 2a through 2d.
 - a. Click the start button on the Windows Taskbar
 - b. Then click the My Computer icon
 - c. Select the D drive under the Devices with removable storage. Double click.
 - d. After that, click the setup icon.
3. Click the Install Holter family of products
4. Under the Welcome to Philips Zymed Holter Family of products installation wizard, click next
5. Review the setup instructions on the Philips InstallShield Wizard and then click next.
6. Review the terms of the License agreement. Select I Accept and click next.
7. Review the information text and click next
8. Confirmation destination folder. Folder should look like the following: C:\Program Files\holter. Select next
9. Only select remote link and documentation. Click next
10. On the current settings section, click next
11. A pop up windows should prompt you to disconnect any Philips recorder connected to the computer. If there is nothing connected, click ok
At this point, install shield will prompt you to insert your hardlock key. Install the key and click ok. (Note: If you are running
12. Windows XP on a newer computer, the install wizard might ask you to run a compatibility test with the software. Click continue anyway)
13. After setup completes, click finish
14. Reboot computer

REMOTE CONFIGURATION SETUP

1. Double click the Remote Link icon on the desktop
2. Click the Tools menu and then select Configure
3. Under the network tab enter your remote code, username, and password provided to you by United Cardiac Services
4. Replace server address with the Central Link Server IP Address 216.70.238.174
5. Under the User Field Labels tab, replace User Field #1 with Institution Name.
6. Under the Mandatory Demographics tab, check the Last Name, First Name, ID Number, Physician, Institution Name, and Scanned by boxes.
7. Click Ok
8. You are now ready to scan
9. If you have ANY problems, call us at **888-310-7374**.